


Manual Nova Comfort


A


	XS	S	M	L	XL
A	69	84	99	110	121
B	46	46	46	80	80
C	80	85	90	95	115

	XS	S	M	L	XL
A	99	114	129	138	148
B	46	46	46	80	84
C	83	88	93	98	115

DA Indhold

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General information

Thank you for choosing a sling from Etac.

Etac's slings is tested by accredited and independent testing institutes and comply with the requirements of the Medical Devices Directive Class 1 products (MDD 93/42/EEC).

The slings also comply with the requirements in standard ISO 10535.

In this manual, the user is the person being lifted.

The assistant is the person operating the sling.

This symbol appears in the manual along with text.

It indicates a need for caution when the safety of the user or assistant may be at risk.

Etac continuously makes improvements to our products. We therefore reserve the right to make changes to our products without prior notice. The dimensions shown in the drawings or other material are for guidance only.

We cannot be held responsible for any printing errors or omissions.

Warranty: Two years for material and manufacturing defects with correct use.

For further information on Etac's range of slings and lifts, see www.etac.com.

Properties

The Nova Comfort sling has a high/low back and undivided leg support.

Approved for users up to 350 kg/770 lbs.

The sling can be used for lifting to and from a lying position or for lifting to and from a sitting position.

Remember

- In order to avoid injury during handling and use, please read through the lift and sling manual.
- It is important to choose the right size.
- It is important to follow the instructions on how to perform a transfer.
- Check the lifting sling regularly, particularly after washing. Pay particular attention to wear and damage to seams, material, straps and loops.
- Do not use damaged lifting accessories.
- Never use a defective lifting sling!
- The sling should not be stored in direct sunlight.
- If required, consider using one or more co-assistants from time to time.
- Plan the lifting operation prior to the lift to ensure it takes place as safely and as smoothly as possible.
- It is important to check that the straps are correctly fastened to the slingbar once the straps are fully stretched, but before the user is lifted from the support.
- Never lift a user higher from the ground than is necessary to perform the lift.
- Never leave a user unattended in a lifting situation.
- The arms are placed on the inside when using a sling with a high back, unlike the low-back slings where the arms can be on the inside or outside.
- If the lifting frame which is used has a catch to prevent accidental unhooking, the catch must be in working condition. Lifting frames without a catch or where the catch does not work must not be used.

EN

Operation

Practical operation

fig. B

Make sure that the wheels of a wheelchair, bed, bunk, etc. are locked during operation.

1. Fit the sling around the user.
Ensure that the arms are in or out depending on the model being used.
2. Suspend the sling from the hooks on the slingbar.

3. Tension the straps of the sling without lifting the user. Check that the sling is correctly suspended from the hooks and is correctly fitted around the user.
4. Lift the user, taking care not to lift him/her any higher than is necessary.

EN

Cleaning

Washable up to 90 °C/194 F.

Tumble-dry up to 45 °C/113 F

NB! Do not bleach, autoclave or dry-clean

The sling must be checked on a regular basis, preferably before each use but especially carefully after washing.

Check there is no wear and damage to seams, material, straps and loops.

NB! Never repair a sling yourself. A worn sling must never be used.

Accessories

Extension strap

Available in lengths 14 cm/
6 Inches or 18 cm/7 Inches, sold
separately in bags.

Marking straps

Marking straps for marking tested
loop positions on the sling for lifting.

Two supplied with every sling in
plastic bags. Can also be purchased
as an accessory.

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Etac Sverige AB
Etac A/S