
DK	FH Strap	2
GB/US	FH Strap	10
DE	Fuß/Hand-Hebeband	20
SE	Fot och hand band	29
FR	Sangle FH	37

Manual – vers. 1.00

CE

Guldmann™

DK FH strap

Vers. 1.00

Varenr.:
28655 FH Strap

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
Vers. 1.00

Item nos:
28655 FH Strap

1.00	Purpose and use	11
1.01	Manufacturer	11
1.02	Purpose	11
1.03	Important/Precautions	11
1.04	Use	11
2.00	Maintenance	13
2.01	Cleaning	13
2.02	The owner's daily maintenance duty	13
2.03	Disposal of sling/accessory	14
3.00	Service and lifetime	14
3.01	Safety/service inspections	14
3.02	Lifetime	16
4.00	Technical specifications	17
5.00	EC-Declaration of conformity	17
6.00	Environmental policy statement – V. Goldmann A/S	17
	WARRANTY – U.S.A. and countries outside the EU	18
A.	Users guide	18
B.	WARRANTY	18
7.00	Fitting the FH Strap	19

1.00 Purpose and use

1.01 Manufacturer

V. Guldmann A/S
Graham Bells Vej 21-23A
DK-8200 Aarhus N
Tel. + 45 8741 3100
Fax + 45 8741 3131
www.guldmann.com

1.02 Purpose

The FH Strap is used for lifting and holding a leg or an arm in hospitals, nursing homes, rehabilitation centres, and in private homes.

The strap is designed for use with both mobile lifters and ceiling hoist systems.

The strap is ideal in connection with operations or when personal hygiene and other patient care tasks are performed.

Conditions for use

- The use of the FH strap is subject to the following:
- The FH strap is used by trained staff or persons who have been instructed in the use of the product in question.
- The FH strap is adjusted in such a way that it encircles the part of the leg or arm to be lifted.
- The maximum nominal load, 100 kg (220 lbs) must not be exceeded.
- The helper pays attention to the well-being of the user when using the FH strap.
- The FH strap is used with the Guldmann lifting hanger

1.03 Important/Precautions

- Read the instructions carefully before using the FH strap.
- The FH strap's maximum load must never be exceeded.
- The FH strap may only be used to lift an arm or a leg .
- Before the FH strap is used, it must be examined according to point 2.02.
- The FH strap is adjusted.in such a way that it encircles the part of the leg or arm to be lifted.
- Possible repairs must only be made by the manufacturer.

1.04 Use

If there is any doubt about the selection or use of a lifting sling/ accessory, please contact your supplier.

Important!


Plan the move. Never leave the leg or arm in the FH strap unattended. Do not start to lift the leg or arm until it has been checked that the user cannot get trapped and that the FH strap does not catch on the bed. Be careful with any tubes and wires that are attached to the user. Check that the hand control and hand control cable is free of hanger, patient and other objects before the hoist is activated up or down moved.

Guldmann shall not be liable for faults or accidents due to incorrect use of the FH strap, or for reasons of inadequate attention on the part of the carer or user. If the FH strap is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.

GH Lifting hanger

Caution!


Be careful when attaching the lifting strap on the hooks. Check that the strap has been pulled completely through the rubber safety catch (A) and into place in the lifting hanger's hooks. When pressing the up button to lift the leg or arm, check again that all the strap remains correctly placed in the lifting hanger's hooks (fig. 1a and fig. 1b).

GH lifting hanger, 4 attachment points

Caution!

Be careful when attaching the lifting strap on the hooks. Check that the strap has been correctly placed in the lifting hanger's hooks. When pressing the up button on the hand control to lift the leg or arm, check again that the strap remains correctly placed in the lifting hanger's hooks (fig. 1).

Placing the FH strap, look at section 7.

Cleaning

Normal washing at the indicated temperature

Do not use bleaching agent

Tumble-drying at low temperature

Do not iron

The owner's daily maintenance duty

Check the lifting sling/accessory for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Sling/accessory inspection checklist

Before using a Guldmann sling/accessory check the following:

Is sling/accessory clean?

Follow facility specific infection control procedure.

Is the sling's/accessory's label present, legible and complete?

Missing, illegible or incomplete label(s) could make identification of appropriate size, function, and or weight limit capacity of the sling/accessory impossible.

Is the fabric intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

Is the fabric intact?

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

Is the sling/accessory the original size and length without the use of knots, pins, tape or other methods to change the shape, shorten or lengthen them?

Conclusion

If the sling/accessory suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

2.03 Disposal of sling/accessory

The slings/accessory are disposed of by incineration. By proper incineration polyester will be degraded to carbon dioxide and water.

3.00 Service and lifetime

3.01 Safety/service inspections

In accordance with international standard EN/ISO 10535 “Hoist for the transfer of disabled persons – Requirements and test methods” an inspection must be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Safe operating practices with sling/accessory

Considerations for damaged or defective sling/accessory and taking them out of service:

Withdraw the sling/accessory from service if one or more of the following conditions are present:

- chemical or caustic burns
- melting or charring of any part of the sling/accessory
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling/accessory tag
- knots in any part of the sling/accessory
- abrasion
- other visible damage that causes doubt as to the strength of the sling/accessory

Sling/accessory inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling/accessory inspection system has additional benefit. Systematic sling/accessory inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling/accessory types and sizes.

Sling/accessory inspection system

Development of a specific procedure and program for the inspection of sling/accessory at your facility is your best safeguard. Consider employing a three part system of inspection. Sling/accessory that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

1) Initial

This level of inspection is done at the time that the sling/accessory is received into your facility. The inspector should ensure that no damage has occurred during transit, and also verify that the sling/accessory work load limits match those contained in the manufacturer's catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage.

2) Frequent

The frequent level of inspection should be done by the sling/accessory user before each use. The sling/accessory should be examined and removed from service if damage is detected. The sling/accessory user should also determine that the sling/accessory is proper for the user conditions, care task required and the required weight capacity.

3) Periodic

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling. Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.

Sling/accessory inspection technique

The sling/accessory inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling/accessory. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting sling/accessory inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the sling/accessory stock number, width and length, the unique sling/accessory identification number (important in differentiating similar slings/accessories), as well as the condition of the products. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

Sample visual examples of synthetic sling/accessory damage ^{x)}


Chemical/caustic burns

Broken stitching


Crushed / Frayed webbing

Knots

Melting / Charring

^{x)} *sample visual images not intended to represent all types of potential damage*

3.02

Lifetime

The life of the product is individual and depends on how it is used, washed etc. Before use the product must be examined according to description in section 2.02 and if it does not meet the inspection requirements, it must be discarded if necessary.

4.00 Technical specifications

Lifting capacity, SWL 100 kg (220 lbs)
Material Polyester

5.00 EC-Declaration of conformity

The product is manufactured in compliance with the Council Directive 93/42/EEC of June 14th 1993, including amendments, as medical device class 1.

6.00 Environmental policy statement – V. Guldmann A/S

Guldmann is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

It is Guldmann's goal to:

- Comply with the current environmental legislation (e.g. WEEE and REACH directives)
- Ensure that we, at the widest possible range, use RoHS compliant materials and components
- Ensure that our products do not have an unnecessary negative impact on the environment regarding use, recirculation or disposal
- Ensure that our products contribute to a positive working environment in the places they are utilised

Inspections are made annually by the Department for Nature and Environment from the Municipality of Aarhus using the Danish Environmental Protection Act, section 42 as a reference.

WARRANTY – U.S.A. and countries outside the EU

A. **Users guide**

Before using the product, read the entire operation manual including warranty.

B. **WARRANTY**

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann designated technician. Any parts or components repaired or replaced by a Guldmann designated technician will be guaranteed for the remainder of the Warranty Period.

The warranty does not cover any part of the equipment which has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment which has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

Service or Repair


Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.


7.00

Fitting the FH Strap


Place the strap on the foot/ankle.

Attach the two lifting straps on the lifting hanger. Lift until all straps are taut, stop and check all attachments are correct.

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Guldmann Srl
Guldmann Inc.