

Flo-Form[®]

Fluidised Positioners



What is Flo-Form?

With its fluidized technology, Flo-Form can be molded to any shape and will hold its form until it is remolded. With its unique characteristic, Flo-Form is an excellent medium for patient positioning and pressure redistribution.



Product Features

- Flo-Form[™] is a fluidized medium, designed to maximize contouring and support
- Effectively positions patients and promotes neutral spine and limb alignment
- The Flo-Form[™] Positioner will not flow on its own, but it effectively contours when normal body forces are applied (regardless of the direction of the force)
- Redistributes pressure over a greater surface area

Suitable for use in many environments



Neonatal



ICU



Elderly Care
Positioning



Disability Care
Positioning


For Head Positioning:

- Redistributes pressure to protect the skin
- Helps reduce lateral skull flattening
- Promotes neutral head/trunk alignment
- Facilitates tube/line management
- Protects the ears, nose and face against shear and friction



For Body Positioning:

- Creates a soft outer boundary to support a flexed, midline position
- Helps maintain neutral body alignment and positive airway
- Encourages normal musculoskeletal and neuromotor development
- Facilitates tube/line management



Flo-Form Fluidized Positioner Range:

Product No.	Product Name
FF3510	Flo-Form Mini - 9x25cm
FF7070	Flo-Form XS Square - 18x18cm
FF7010	Flo-Form X-Small - 18x25cm
FF9015	Flo-Form Small - 23x38cm
FF1220	Flo-Form Medium - 30x50cm
FF1824	Flo-Form Large - 45x60cm
FF2028	Flo-Form X-Large - 50x70cm
FT0730	Flo-Form Tube 18x75cm
FT0740	Flo-Form Tube 18x100cm
FT0755	Flo-Form Tube 18x140cm
FF0400	Flo-Form Medium with blue cover - 30x45cm
FF0100	Flo-Form Large with blue cover - 50x64cm



- Clean topically using hospital-approved products
- DEHP and latex free