
Auf der Gröb 2 c 83064 Raubling Germany

Tel.: +49 (0)8035-9580448 Fax: +49 (0)8035-9580446

info@hernik.de • www.hernik.de





	1	0	rd	er
_	_	$\overline{}$		

Offer

Helmet CAPOVARIO

Exportprices Release: 2023

BASIC

technical specification:

- solid, anatomic frame made of synthetical, size adjustable
- molded rigid foam in 3D-concavity
- microfiber: washable by 40 °C, can be disinfected, breathable, flame resistant, antistatic, dirt- and water-repellent
- forehead protection included
- upholstery inlay (microfiber),
- weight 280 320 g
- all accessories can be adapted later, too!







(60055Y-Z)

size (Y)	0	1	2	3	4
circumference of the head	44 - 48 cm	48 - 51 cm	52 - 56 cm	56 - 60 cm	58 - 62 cm
weight	270 g	280 g	290 g	305 g	320 g
PINKRED	☐ 6005X0-P	☐ 6005X1-P			
TURQUOISE BLUE	☐ 6005X0-T	☐ 6005X1-T			
RED			☐ 6005X2-R	☐ 6005X3-R	☐ 6005X4-R
BLACK			☐ 6005X2-S	☐ 6005X3-S	☐ 6005X4-S
GREY			☐ 6005X2-G	☐ 6005X3-G	☐ 6005X4-G
BLUE			☐ 6005X2-B	☐ 6005X3-B	☐ 6005X4-B
Strap Please choose ONE of them!	□ standard (600530)		padded, for size 0, 1, 2 (600535-1)padded, for size 3, 4 (600535-2)		

colours:













PINKRED RED BLACK GREY BLUE



	accessories						
☐ 4 strapes (600590-X) ☐ size 0 ☐ size 1 ☐ size 2 ☐ size 3 ☐ size 4	closed top side soft	FIDLOCK One-Hand-Use (600585-X)					
upholstery for the buckle (600580-X) 1: for standard strap 2: for padded strap protects the skin - no contact with the buckle	magnetic lock (600550-X) 1: for standard strap 2: for padded strap note: • user is not able to put away the helmet • a special magnetic key must be used!	chin protector (600515-X) 1: for size 0/1/2 2: for size 3/4					
ear protectors* (600510-X) 1: for size 0/1/2 2: for size 3/4 *padded strap necessary!	visor# (600570-X) 1: for size 0/1/2 2: for size 3/4 # delivery time 2 - 3 weeks	cheek protector (600565-X) 1: for size 0/1/2 2: for size 3/4					
spare parts							
additional strap, padded (600535-X) 1: for size 0/1/2 2: for size 3/4	alternating inlay, 2 straps (600540-X) size 0 size 1 size 2 size 3 size 4	alternating inlay, 4 straps (600585-X) size 0 size 1 size 2 size 3 size 4					
production in two steps (600545-X) A blanks can be adjusted with hot air by a technician to adapt it on place to the individual requirements of the patient's individual requirements before finally obtained and finished in our workshop.	Fixation of all movable elements (600555-X) After a successful fitting on place, all removable parts are fixed in our workshop. This may be important for very vigorous patients.	mouth protection, adjustable* (600560-X) 1: for size 1/2 2: for size 3/4					

Washing then only by hand!

*padded strap necessary!