



## Hand pump

The vacuum hand pump is used to form small and medium cushions. Simple and reliable, it comes included with each small-medium cushion.



## Electro Vacuumer™ - Electric pump

The vacuum electric pump makes it easier to shape large Stabulo forms thanks to its high performance. The 12V power supply ensures safety of use and enables connection to the car cigarette lighter socket. It comes with Grande mattresses and Beanseat.



## Neoprene strap

The elastic neoprene strap helps to better stabilize the user in a cushion. Thanks to the Velcro surface, it can be fastened to the cushion at various heights. Each Comfortable Plus and Comfortable Duo model is equipped with such a strap as standard.



## Mounting strap

Universal, adjustable stroller strap is used to mount the cushion in a wheelchair, armchair or chair. It prevents it from moving and thus makes it easier to place the patient into the cushion.



## Shorts stabilizing the pelvis and provides abduction of thighs

Shorts are used to stabilize the pelvis and provides abduction of thighs. They can be used to position the pelvis and help restore proper posture. They offer freedom of upward movement of hips, while blocking forward movement of pelvis, which is advisable in case of large spastic conditions. With shorts you can better position your hips in Beanseat of Velcro type.

