**QS7X**

OPTIONS AND ACCESSORIES CATALOGUE



Errors and omissions excepted - Information in this product brochure may contain unintentional errors or inaccuracies. All weights shown in this catalog are for reference only. They are calculated using Seat Width (SW) 400 mm/Seat Depth (SD) 420 mm dimensions. Please note that in practice the actual weight may vary from the theoretical weight due to variations in manufacturing tolerances and SW/SD dimensions. The information contained herein is based on our present knowledge and experience and is given in good faith. However, variations in weight are not grounds for a refund.

Seat Slings

Nylon, tension adjustable

Item Code: QSF020153, QSS050126

Nylon seat sling with adjustable tension straps for individualized support and seating stability.

309g



Seat Slings

1 utility bag

Item Code: QSF020154, QSS020154

One utility bag underneath the seat sling - useful for storing belongings and compact tools.

41g




Seat Slings

2 utility bags

Item Code: QSF020155, QSS020155

Two utility bags underneath the seat sling - useful for storing belongings and compact tools.

82g




Seat cushion, cover with zipper, black

Foam (medium-soft) Thickness: 30 mm

Item Code: QSF020101, QSS020101

Simple foam cushion, available in 30 mm height (020101) and 50 mm height (020102).

625g




Seat cushion, cover with zipper, black

Foam (medium-soft) Thickness: 50 mm

Item Code: QSF020102, QSS020102

Simple foam cushion, available in 30 mm height (020101) and 50 mm height (020102).

870g



JAY Cushion

JAY Basic

Item Code: JAY020002

Mild contoured foam cushion, curved foam base for use with sling seat.

549,5g




JAY Cushion

JAY Soft Combi P

Item Code: JAY020003

Pre-contoured foam base with lateral and medial thigh support, incontinence resistant coated, easy-clean.

640,5g



JAY Cushion

JAY Easy Visco C

Item Code: JAY020001

Pre-contoured foam base with medial/lateral thigh support and a Visco elastic foam in the seat well for immersion of ischial tuberosities choice of curved or flat base.

985,5g



JAY Cushion

JAY Easy Fluid

Item Code: JAY020007

Pre-contoured foam base, built-in medial and lateral thigh support, choice of curved or flat base, JAY Flow fluid Tripad.

1750g



JAY Cushion

JAY LitePro

Item Code: JAY020067

One of the lightest, low-profile adjustable skin protection cushions.

Designed for active self-propellers, it offers adjustable pressure management, improved breathability, and stable support.

920g



Seat accessories

Lap belt with buckle, metal

Item Code: QSF090018, QSS090018

Pelvic positioning belt, that can prevent sliding and falling out of the wheelchair. (Pair of weight)

250g



Push handles

Push handles, long

Item Code: QSF030201, QSS030201

Standard handles provide optimal control when an attendant is pushing the chair but still keep the chair's overall length compact. (Pair of weight)

384g



Push handles

Push handles, fold-down

Item Code: QSF030203, QSS030203

The fold-down push handles may be used by pressing the button inside. Although they have less ease of grip than standard push handles, they are lightweight and have a discreet profile that nearly disappears when folded. (Pair of weight)

462g



Push handles

Push handles, height-adjustable

Item Code: QSF030204, QSS030204

Height adjustable push handles allow adjustment to a comfortable height that can be set independent of the backrest height. They can be dropped down and tucked away when not in use or be fully removed if they are no longer needed. (Pair of weight)

850g



Push handles

Push handles, height-adjustable for JAY backs

Item Code: QSF030205, QSS030205

Height adjustable push handles for JAY backrest. Use these if the chair also has a JAY backrest installed. These push handles can be utilized in combination with a low back height or short backrest to achieve desired push handle height for caregivers without interfering with the users propulsion. They can be dropped down and tucked away when not in use or be fully removed if they are no longer needed. (Pair of weight)

850g



Without push handles (L)

Item Code: QSF030200, QSS030200

Wheelchair comes without push handles.

162g




Backrest, fixed (90°)

Item Code: QSF030010, QSS030010

Most rigid backrest option due to fewer moving parts. This option should only be selected if the user is established and unlikely to change their upholstery. This option does not have angle adjustment for fine tuning or changing needs. Backrest tubes are perpendicular to the seat. Weight is included in the basic chair weight.

-50g




Desk sideguard preparation for later replacement

Item Code: QSF030029, QSS030029

Uses a specific backrest tube that allows the desk sideguard to be retro-fitted at a later stage. Without this option, the backrest tube would need to be replaced to enable installation.

44g



Backrest, fold-down (to the front)

Item Code: QSF030031, QSS030031

Backrests that fold forward are designed to reduce the height of your wheelchair to fit into everyday cars and to reduce the overall package size for stowage.

1080g



Backrest, angle-adjustable

Item Code: QSF030012, QSS030012

Back angle is adjustable forwards and backwards, (range is dependant on the chair (see order form)) giving the ability to squeeze the client when a deep bucket is used or open the backrest for those with limited hip ROM.

650g



Backrest upholstery

EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black

Item Code: QSF030316, QSS030316

The EXO backrest is designed for both comfort and breathability. The ribbed curves allow for air to flow and the compressed EVA fabric ensures that it is both lightweight and comfortable.

547g




Backrest upholstery

EXO PRO Backrest sling, padded, adjustable, breathable, spacer fabric, 1 accessory bag, black

Item Code: QSF030317, QSS030317

The EXO PRO combines the features of the EXO upholstery with added 3DX spacer fabric technology providing additional comfort and extreme breathability.

547g



Backrest accessories

Stabilising bar, swing-away

Item Code: QSS030401

Swing-away stabilising bar designed to enhance frame rigidity during use, while allowing it to be moved aside for easier transfers and access.

338g



Backrest accessories

Headrest, height-, depth- and angle-adjustable medium (200 x 110 mm)

Item Code: QSS030411

A headrest provides support and stability for the head and neck for individuals who have limited head and neck control. Available in tall height (030410) and medium height (030411)

1350g



Backrest accessories

**Headrest, height-, depth- and angle-adjustable tall
(230 x 110 mm)**

Item Code: QSS030410

A headrest provides support and stability for the head and neck for individuals who have limited head and neck control. Available in tall height (030410) and medium height (030411)
1390g



Backrest accessories

Stabilising bar, automatic folding

Item Code: QSF030402, QSS030402

The stabiliser bar provides extra rigidity to wider wheelchairs with upholstery rather than solid backrests. Folds up automatically when the chair is folded.
260g



JAY J3 Back - Shallow Contour (S) Upper Thoracic - 420 mm

Item Code: JAY030012

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.
2400g



JAY J3 Back - Shallow Contour (M) Mid Thoracic - 340 mm

Item Code: JAY030010

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.
1670g



JAY J3 Back - Shallow Contour

(M) Upper Thoracic - 460 mm

Item Code: JAY030019

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

2520g



JAY J3 Back - Shallow Contour

(T) Mid Thoracic - 380 mm

Item Code: JAY030033

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

1830g



JAY J3 Back - Shallow Contour

(T) Upper Thoracic - 500 mm

Item Code: JAY030025

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation. The image shows option JAY030042 (460 mm).

2620g



JAY J3 Back - Mid Contour

(S) Upper Thoracic - 420 mm

Item Code: JAY030013

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

2460g



JAY J3 Back - Mid Contour

(M) Mid Thoracic - 340 mm

Item Code: JAY030011

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

2240g



JAY J3 Back - Mid Contour

(M) Upper Thoracic - 460 mm

Item Code: JAY030020

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

2640g



JAY J3 Back - Mid Contour

(T) Mid Thoracic - 380 mm

Item Code: JAY030018

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

2340g



JAY J3 Back - Mid Contour

(T) Upper Thoracic - 500 mm

Item Code: JAY030026

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

2680g



JAY J3 Back - Mid Deep Contour (S) Upper Thoracic - 420 mm

Item Code: JAY030045

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

3120g



JAY J3 Back - Mid Deep Contour (M) Mid Thoracic - 340 mm

Item Code: JAY030044

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

1650g



JAY J3 Back - Mid Deep Contour (M) Upper Thoracic - 460 mm

Item Code: JAY030042

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

3250g



JAY J3 Back - Mid Deep Contour (T) Mid Thoracic - 380 mm

Item Code: JAY030041

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

1720g



JAY J3 Back - Mid Deep Contour (T) Upper Thoracic - 500 mm

Item Code: JAY030046

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

3380g



JAY J3 Carbon Back - Shallow Contour

Short (S) Mid Thoracic - 170 mm

Item Code: JAY060001

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

980g



JAY J3 Carbon Back - Shallow Contour

Short (S) Upper Thoracic - 300 mm

Item Code: JAY060003

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

980g



JAY J3 Carbon Back - Shallow Contour

Tall (T) Mid Thoracic - 240 mm

Item Code: JAY060002

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

1040g



JAY J3 Carbon Back - Shallow Contour

Tall (T) Upper Thoracic - 380 mm

Item Code: JAY060004

The JAY J3 wheelchair back product line offers a multitude of width, height, and contour depth back shells to fit the back to the rider. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

1240g



Whitmyer Headrest

COBRA XTRA, flip back

Item Code: WHT030415

Cobra Xtra Flip Back Hardware - sleek, curved, and modular design with superior positioning capabilities, sizes ranging from early intervention to adult, and flip back function. Plush Pad - provides posterior support with a high degree of comfort and safety. This headrest is ideal for tilt or recline, accepts additional lateral and anterior supports, and is easily adjustable to adapt with clients as their needs change.

1534g



Whitmyer Headrest

AXYS, Standard

Item Code: WHT030416

Axys hardware - durable, compact, one-size option with basic anterior- posterior adjustability at an outstanding value. Plush Pad - provides posterior support with a high degree of comfort and safety. This headrest is ideal for tilt or recline, accepts additional lateral and anterior supports, and is easily adjustable to adapt with clients as their needs change.

1522,5g



Sideguard

Carbon fibre, cold weather protection, slim (black)

Item Code: QSF040101, QSS040101

Ultra light weight and strong carbon sideguards protect against dirt and moisture and keep the cushion in place. Inner side has narrow foam for padding and cold weather protection. They are non removable. (Pair of weight)

416g



Sideguard

Aluminium, cold weather protection (in frame colour)

Item Code: QSF040102, QSS040102

Aluminium sideguards protect against dirt and moisture and keep the cushion in place. Inner side has narrow foam for padding and cold weather protection. They are non removable. (Pair of weight)
780g



Sideguard

Aluminium, fender, cold weather protection (in frame colour)

Item Code: QSF040103, QSS040103

Aluminium sideguards with fender provide the best protection against dirt and moisture and keep the cushion in place. Inner side has narrow foam for padding and cold weather protection. They are non removable. (Pair of weight)
910g



Sideguard

Carbon fibre EVO, fender, cold weather protection (black)

Item Code: QSF040104, QSS040104

Ultra light weight and strong carbon sideguards protect against dirt and moisture and keep the cushion in place. Inner side has narrow foam for padding and cold weather protection. They are non removable. (Pair of weight)
510g



Sideguard

Composite, removable, 210 mm high (310 mm long), black

Item Code: QSF040107, QSS040107

Composite sideguards protect against dirt and moisture and keep the cushion in place. This style is detachable for side transfers and lowering transport weight. (Pair of weight)
890g



Sideguard

Composite, fender (black)

Item Code: QSF040120, QSS040120

Composite sideguards protect against dirt and moisture and keep the cushion in place. This style is non detachable. (Pair of weight)
966g



Desk sideguard

Flip-up, removable

Item Code: QSF040005, QSS040005

This armrest offers arm support with a built in side guard protection and flips away to allow transfers into the wheelchair. (Pair of weight)
1828g



Desk sideguard

Armpad (250 mm), fixed, flip-up, removable

Item Code: QSF040001, QSS040001

This armrest offers arm support with a built in side guard protection, arm pad lengths of 250 mm, the height of this armpad cannot be altered. (Pair of weight)
600g



Desk sideguard

Armpad (250 mm), height-adjustable, flip-up, removable

Item Code: QSF040002, QSS040002

This armrest offers arm support with a built in side guard protection, arm pad lengths of 250 mm, the height of this armpad can be altered. (Pair of weight)
700g



Desk sideguard

Armpad (310 mm), fixed, flip-up, removable

Item Code: QSF040003, QSS040003

This armrest offers arm support with a built in side guard protection, arm pad lengths of 310 mm, the height of this armpad cannot be altered. (Pair of weight)

750g



Desk sideguard

Armpad (310 mm), height-adjustable, flip-up, removable

Item Code: QSF040004, QSS040004

This armrest offers arm support with a built in side guard protection, arm pad lengths of 250 mm, the height of this armpad can be altered. (Pair of weight)

850g



Single-post sideguard

Armpad (330 mm), depth adjustable, height-adjustable, removable

Item Code: QSF040052, QSS040052

Our lightest armrest offers the most support and protection with a built-in sideguard. It can be adjusted and moved forward/ backward, as well as being height adjustable. (Pair of weight)

1934g



Tubular armrest. Alu, swing away, padded, removable

Item Code: QSF040117, QSS040117

Armrest (composite side guards can be ordered separately). This style is detachable for side transfers and lowering transport weight. (Pair of weight)

1522g



Hanger

Aluminium, swing-away, LH & RH

Item Code: QSS05006, QSS05007

High performance aluminum hanger with swing-away mechanism for easier transfer, lightweight. Available in left (050006) and right (050007).
270g



Tray

Push-fit, SW 380 - 480 mm

Item Code: QSS090011

Tray attachable to the wheelchairs, providing a firm table for eating, working, reading and other activities. Push fit mechanism.
2320g



Hanger

Hanger angle 80° (100°)

Item Code: QSS050001

80° hanger angle positions feet closer to wheelchair frame for a tighter turning circle and more support directly through the feet. (Pair of weight, aluminium hanger on image)
595g



Hanger

Hanger angle 70° (110°)

Item Code: QSS050002

70° hanger angle for a more extended knee position where this is required for comfort or postural needs. (Pair of weight, aluminium hanger on image)
605g



Hanger

Hanger angle 90°

Item Code: QSS050041

90° hanger angle positions feet closer to wheelchair frame for tightest turning circle and more support directly through the feet. (Pair of weight) 585g



Hanger

Hanger, aluminium, flip-up (0° - 90°), swing-away left & right

Item Code: QSS050009 - QSS050010

Flip-up aluminum hanger with rigid aluminum housing. Easy operation, swing-away mechanism for easier transfer and durability (minimised wear & tear). Available with left-mount (050009) and right-mount (050010).

1179g



Footboard

Divided, composite, angle-adjustable, flip-up (side-ways), heel loop, black

Item Code: QSS050020

Divided composite footboard, that is both angle and depth adjustable, and flips up sideways for easy transfer. Includes heel loop that can be utilized to assist with lower extremity positioning. (Pair of weight) 656g



Footboard

Divided, aluminium, angle-adjustable, flip-up (side-ways), heel loop, black

Item Code: QSS050050

Divided footboard, made of aluminum, angle and depth adjustable, flip up sideways for easy transfer. (Pair of weight) 1820g



Footboard

Platform, aluminium, angle-adjustable, flip-up (side-ways), heel loop, black

Item Code: QSF050102, QSS050102

Aluminum platform footboard, angle-adjustable for better positioning, flip-up (sideways) for easier transport, heel loop keeps foot in place.
858g



Swing-away mechanism for platform footboard

Item Code: QSF050126, QSS050126

Swing-away mechanism for Alu platform footboard. The weight of this option is added to option XXX050102. The option is only available in combination with it.
45g



Locking system, platform footplate, aluminium

Item Code: QSF050062, QSS050062

Locking mechanism in the middle for increased stability and to prevent unintended opening of the footboard. Unlocks for folding. The weight of this option is added to option XXX050102. The option is only available in combination with it.
36g



Footboard

Platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap

Item Code: QSF050059, QSS050059

Lightweight aluminium platform footboard with adjustable angle setting and sideways flip-up function for easier transfers. Includes calf strap for secure lower leg positioning.
598g



Footboard

Platform, lightweight, plastic, angle & depth-adj. flip up (sideways), calf strap, black

Item Code: QSF050131, QSS050131

Composite footboard with angle and depth adjustment for full foot support.

716g



Footboard

Platform, lightweight, carbon fibre, angle & depth-adj. flip up (sideways), calf strap, black

Item Code: QSF050132, QSS050132

Carbon fibre footboard that flips up sideways with angle and depth adjustment for full foot support.

676g



Footboard

Platform, lightw., plastic, angle & depth-adj. flip up (sideways), high-mount

Item Code: QSF050061

Lightweight plastic platform footboard with adjustable angle and depth settings. Sideways flip-up design and high-mount option support compact setup and easier transfers.

908g



Footboard

Platform, lightweight, carbon fibre, auto-folding, angle & depth-adjustable, calf strap, black

Item Code: QSF050027

Carbon fibre footboard with angle and depth adjustment for full foot support. Auto-folds for easy folding.

715g



Adjustable side protection, to create a side-stop on the footboard

Item Code: QSF050127, QSS050127

Lateral protection reduces risk of feet sliding sideways from the foot support.

40g



Frame 100 °

Item Code: QSF010012

When selecting frame angle, lower extremity position is crucial. This depends on the angle of the lower leg support (e.g., hanger or front rigging) to support thighs and lower legs at the desired angle. Frame angle impacts length, turning radius, and the user's ability to approach objects. (Weight included in basic chair weight.)

0g



Frame 92 ° , Performance

Item Code: QSF010013

When selecting frame angle, lower extremity position is crucial. This depends on the angle of the lower leg support (e.g., hanger or front rigging) to support thighs and lower legs at the desired angle. Frame angle impacts length, turning radius, and the user's ability to approach objects. (Weight included in basic chair weight.)

0g



Frame 88 ° , Performance Compact

Item Code: QSF010009

When selecting frame angle, lower extremity position is crucial. This depends on the angle of the lower leg support (e.g., hanger or front rigging) to support thighs and lower legs at the desired angle. Frame angle impacts length, turning radius, and the user's ability to approach objects. (Weight included in basic chair weight.)

0g



Frame versions

Open frame - Basic frame chair weight

Item Code: QSF010067

The open fixed front frame simplifies car transfers. Fully adjustable for a personalised setup. Basic chair weight - Base frame of chair including backrest tubing (tubes without PH QSF030200, QSS030200 Without push handles (L)), castor links, & seat upholstery 4538g



Frame versions

Pro - Difference from Open frame basic frame chair weight

Item Code: QSF010068

The Pro frame features a fixed welded backrest and welded axle stems, reducing moving parts, cutting additional 300 g and maximises rigidity. Recommended for experienced riders who do not require later adjustment. Base frame of chair including backrest tubing (tubes without PH QSF030200, QSS030200 Without push handles (L)), castor links, & seat upholstery 4098g (-440g)



Frame versions

Hybrid - Difference from Open frame basic frame chair weight

Item Code: QSF010035

The reinforced fixed front frame delivers exceptional rigidity and supports users up to 140 kg. Its strength makes it ideal for use with power assist devices, while steeply angled tubes allow independent car transfers. Base frame of chair including backrest tubing (tubes without PH QSF030200, QSS030200 Without push handles (L)), castor links, & seat upholstery 5007g (+469g)



Frame versions

Folding front frame - Difference from Open frame basic frame chair weight

Item Code: QSF010070

With the new Front-Fold mechanism, the front frame, caster link, and fork fold in one smooth motion minimizing QS7 X's size for an exceptionally compact package that's easy to store and transport. Base frame of chair including backrest tubing (tubes without PH QSF030200, QSS030200 Without push handles (L)), castor links, & seat upholstery 4910g (+372g)



Frame versions

Swing Away - Difference from Open frame basic frame chair weight

Item Code: QSS010015


Features footrests that move to the side, making transfers more comfortable, especially from standing. Removing the footrests lets you position the wheelchair closer for easier transfers and reduces the overall pack size. Base frame of chair including backrest tubing (tubes without PH QSF030200, QSS030200 Without push handles (L)), castor links, & seat upholstery
4990g (+452g)



Inset: Standard (0 mm)

Item Code: QSF010100


Front frame with 0 mm inset, providing a neutral lower leg alignment and balanced overall setup. Supports a natural knee position and standard wheelchair length.
0g



Inset: Slight (10 mm)

Item Code: QSF010101


Front frame with 10 mm inset, bringing the caster wheels closer to the frame. This results in improved maneuverability, a more compact footprint, and more efficient propulsion energy, while supporting balanced lower leg alignment.
0g



Twin-LED light (mounted on caster stems, battery and controller on calf-strap, details on backpage)

Item Code: QSF110156, QSS110156

Twin-LED light increases visibility and ensures that paths and roads are easier to see even in the dark. Increases visibility at night. (Pair of weight)
SA= 90g ; FF= 212g



Design Frame Protection

Item Code: QSF010040

Protective elements designed to help reduce wear and damage to the wheelchair frame during everyday use and transfers.

40g



Transport Protection, Black

Item Code: QSF010069

Designed to help reduce the risk of damage to the wheelchair during transport and handling.

Lightweight 30g, Standard 46g



Axle plate Adapter, standard

Item Code: QSF070703, QSS070703

Standard axle plate adapter allowing camber adjustment of 0°, 2°, or 4°. Enables quick center of gravity adjustment (± 1 position, approx. 10 mm depending on seat angle). Larger and heavier than the lightweight adapter.

484g



Axle plate Adapter, lightweight

Item Code: QSF070701, QSS070701

Compact and lightweight axle plate adapter with fixed 0° camber. Shorter and smaller design for reduced overall weight.

224g



Axle plate, aluminium, long (full adjustment)

Item Code: QSF010411, QSS010411

Long axle plate version offering full adjustability of the rear seat height for precise setup and positioning.

568g



Axle plate, aluminium, shorted (max. adjustment rear seat height +10mm)

Item Code: QSF010412, QSS010412

Shortened axle plate version, trimmed to allow a maximum rear seat height increase of 10 mm. Reduces weight and offers a cleaner appearance - especially beneficial for low rear seat height configurations.

Maximum weight saving: 166 g (26 g per step)



Frame Accessories

Tie down brackets for transportation according to ISO 7176-19 (max. user-weight 136 kg)

Item Code: QSF090100, QSS090100

Tie-down brackets compliant with ISO 7176-19 for secure wheelchair transportation in vehicles. Approved for a maximum user weight of 136 kg.

292g



Frame Accessories

Crutch holder, loop

Item Code: QSF090009, QSS090009

Loop-style holder designed to securely position crutches on the wheelchair for convenient transport.

125g



Frame Accessories

Anti-tip tubes, plug-in, pair

Item Code: QSF090050, QSS090050

Anti-tip tubes increase wheelchair safety by reducing the risk of rear tipping.

644g



Frame Accessories

Tip assist

Item Code: QSF090001-QSF090003, QSS090001-QSS090003

For carers to overcome kerbs easily.

215g



Frame Accessories

Anti-tip tubes, swing-away

Item Code: QSF090004 - QSF090006 , QSS090004 - QSS090006

Anti-tip tubes increase wheelchair safety by reducing the risk of rear tipping. Available in left (090004), right (090005) and pair (090050).

1084g



Frame Accessories

Transit wheels

Item Code: QSF090010, QSS090010

Transit wheels allow removal of rear wheels for short periods of time to access very narrow spaces, such as a very small doorway or aisle. Please note, the weight must be fully off the rear wheels for removal. (Pair of weight)

920g



Bundling

With F55 14" order, not mounted to the wheelchair

Item Code: EM0000101

QS7 X comes together with Empulse F55. F55 not mounted on the wheelchair.

13800g



Bundling

With F35 order, not mounted to the wheelchair

Item Code: EM0000301

QS7 X comes together with Empulse F35. F35 not mounted on the wheelchair.

10028g



Bundling

With R90 order, not mounted to the wheelchair

Item Code: EM0000401

QS7 X comes together with Empulse R90. R90 not mounted on the wheelchair.

8900g



Bundling

With M90 order, not mounted to the wheelchair

Item Code: EM0000501

QS7 X comes together with Empulse M90. M90 not mounted on the wheelchair.

13992g



Castor fork

Carbotecture® , (105 mm) Colour: Black

Item Code: QSF080013, QSS080013

Lightweight 105 mm castor fork made from Carbotecture® material in black, combining reduced weight with high strength and durability.
163g



Castor fork

Aluminium, Black (115 mm) (anodized)

Item Code: QSF080029, QSS080029

Castor fork, aluminium, black (anodized) (Pair of weight)
115mm = 317g ; 90mm = 134g ; 131mm = 350g ; 166mm = 450g



Castor fork

Aluminium, Powder coated (115 mm)

Item Code: QSF080064, QSS080064

Castor fork aluminium powder coated (Pair of weight)
115mm = 317g ; 90mm = 134g ; 131mm = 350g ; 166mm = 450g



Castor fork

Aluminium, one-arm (105 mm) Colour: Black

Item Code: QSF080008, QSS080008

Single arm forks decrease the front frame footprint slightly and adds an element of style to the front of the chair. (Pair of weight)
300g




Castor fork

Aluminium, suspension (Frog Legs ®), black

Item Code: QSF080007, QSS080007

Aluminium castor fork with integrated suspension system from Frog Legs, designed to help absorb vibrations and improve ride comfort.

580g



Castor

Pneumatic, rim: composite, black

Item Code: QSS080305

Pneumatic tires compensate for uneven ground better and offer a more cushioned ride. (Pair of weight)

7" = 868g



Castor

Solid (hard), rim: composite, black

Item Code: QSF080301, QSS080301

Firm castor wheels have lower rolling resistance indoors. (Pair of weight)

4" = 502g ; 6" = 660g ; 7" = 1204g



Castor

Solid (soft), rim: composite, black

Item Code: QSF080302, QSS080302

Width and softer material of these castor wheels provides better grip and rolling over pavements and other everyday outdoor surfaces but slightly more rolling resistance on smooth, flat surfaces. (Pair of weight)

5" = 823g ; 6" = 1138g



Castor

Solid (hard), light-up, rim: composite, black

Item Code: QSF080303, QSS080303

Casters with integrated light elements, lower rolling resistance compared to soft solid casters, may be less comfortable on uneven terrain than soft solid casters. (Pair of weight)

4" = 338g ; 5" = 422g



Castor

Solid (hard), rim: aluminium Colour, rim: Silver

Item Code: QSF080304, QSS080304

Lower rolling resistance compared to soft solid casters, may be less comfortable on uneven terrain than soft solid casters. (Pair of weight)

4" = 343g ; 5" = 442g



Castor

Solid (soft), rim: aluminium Colour, rim: Silver

Item Code: QSF080306, QSS080306

Width and softer Poly-Urethan (PU) of these castor wheels provides better grip and rolling over pavements and other everyday outdoor surfaces but slightly more rolling resistance on smooth, flat surfaces. (Pair of weight)

4" = 428g ; 5" = 602g



Castor

Solid (soft), rim: carbon-composite, ultralight

Item Code: QSF080060, QSS080060

Ultralight solid castor with carbon-composite rim and softer Polyurethane (PU) tread. The increased width and softer PU provide improved grip and smoother rolling over pavements and everyday outdoor surfaces, with slightly higher rolling resistance on smooth, flat floors.

4" = 312g



Rear wheel

Universal, 36 spokes (silver), cross-wise spokes

Item Code: QSF070001, QSS070001

Robust, low-maintenance standard wheels suitable for heavy use and a variety of environments. (Pair of weight, weight specification for 24 inch) 1670g



Rear wheel

Design, 36 spokes (silver), straight spokes

Item Code: QSF070002, QSS070002

Appealing design with stylish, highly polished aluminium, radial spokes and low-profile wheel-hub. (Pair of weight, weight specification for 24 inch) 1586g



Rear wheel

Lightweight, 24 spokes (silver), straight spokes

Item Code: QSF070003, QSS070003

Appealing design with a reduced number of spokes for a discreet appearance. (Pair of weight, weight specification for 24 inch) 1222g



Rear wheel

Proton, 24 spokes (black), straight spokes

Item Code: QSF070004, QSS070004

These are our lightest wheels and feature a sleek look with black wheels and hub and proton branding. Due to their featherweight construction they have a lower user weight rating than other options (100kg weight-limit). (Pair of weight, weight specification for 24 inch) 1150g



Rear wheel

Spinergy, 18 spokes (black), straight spokes; hub: silver

Item Code: QSF070005, QSS070005

Made with a double-walled custom design aluminium rim for strength. Very little maintenance required. Available with 18 spokes (070005) or 12 spokes (070007). (Pair of weight, weight specification for 24 inch) 1390g



Rear wheel

Spinergy, 12 spokes (black), straight spokes, hub: black

Item Code: QSF070007, QSS070007

Made with a double-walled custom design aluminium rim for strength. Very little maintenance required. Available with 18 spokes (070005) or 12 spokes (070007). (Pair of weight, weight specification for 24 inch) 1320g



Rear wheel

Drum brakes, 36 spokes (silver), cross-wise spokes

Item Code: QSF070008, QSS070008

This wheel option is compatible with our drum brake accessory. (Pair of weight, weight specification for 24 inch) 2020g



Rear wheel

Mountainbike, 36 spokes (silver), straight spokes

Item Code: QSF070009, QSS070009

This wheel comes complete with a mountain bike tyre ready for off road use. (Pair of weight, weight specification for 24 inch) 2420g



Quick-release axle

Item Code: QSF070150, QSS070150

Allows rear wheels to be removed quickly and easily for reduced transport size and simplified handling.

224g



Quick-release axle superlight, Titan

Item Code: QSF070153, QSS070153

Lightweight titanium quick-release axle for reduced overall weight and easy rear wheel removal during transport and storage.

134g



Quick-release axle, designed for quads

Item Code: QSF070151, QSS070151

Specially designed quick-release axle with extended grip surface to support users with limited hand function, enabling easier and more secure rear wheel removal.

234g




Tyre

Pneumatic, fine profile, Schwalbe Right Run

Item Code: QSF070101, QSS070101

Classic fully pneumatic tyre for a cushioned ride and minimal rolling resistance. (Pair of weight, weight specification for 24 inch)

800g



Tyre

Pneumatic, slick, Schwalbe One

Item Code: QSF070107, QSS070107

This tyre features slick tread for manoeuvrability and lightweight pneumatic construction with a Kevlar liner to resist punctures. (Pair of weight, weight specification for 24 inch)
600g



Tyre

Puncture proof, tread (+)

Item Code: QSF070102, QSS070102

The free solid tyre option provides the ultimate maintenance-free durability over the life of the tyre. Although these tyres won't go flat, they offer less cushioning and performance compared to pneumatic options. (Pair of weight, weight specification for 24 inch)
1150g



Tyre

Pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution

Item Code: QSF070104, QSS070104

The extremely popular Marathon Plus Evolution provides a practical compromise between the flat prevention of a solid and the light, cushioned performance of a pneumatic tyre. It features a thick solid foam liner with an inner pneumatic section. (Pair of weight, weight specification for 24 inch)
1240g



Handrims

Aluminium, powder-coated, clear

Item Code: QSF070201, QSS070201

Classic lightweight aluminium handrim. Powder coating adds grip especially for users with impaired dexterity or hand function. (Pair of weight, weight specification for 24 inch)
690g



Handrims

Aluminium, anodised, silver

Item Code: QSF070202, QSS070202

Classic lightweight aluminium handrim slides easily under the users hands and resists the appearance of scratches. (Pair of weight, weight specification for 24 inch)

610g



Handrims

Aluminium, powder-coated, coloured

Item Code: QSF070203, QSS070203

Classic lightweight aluminium handrim with the option to colour match to your favourite colour. (Pair of weight, weight specification for 24 inch)

690g



Handrims

Stainless steel, silver

Item Code: QSF070204, QSS070204

Durable stainless steel option that resists the appearance of scratches but is heavier than other handrim options. (Pair of weight, weight specification for 24 inch)

1370g



Handrims

Titanium, untreated

Item Code: QSF070205, QSS070205

Durable titanium handrims with untreated surface, offering high strength, low weight, and a direct grip feel for efficient propulsion.

803g



Handrims

Titanium, powder-coated, clear

Item Code: QSF070206, QSS070206

Lightweight titanium handrims with clear powder-coated finish, combining high strength with added surface protection and a smooth grip feel. 823g



Handrims

Supergrip®

Item Code: QSF070207, QSS070207

This handrim has a rubberised grip that may be helpful for people with upper limb weakness or limited hand function, since it does not slide as easily as other handrim options when a user is slowing the wheelchair. (Pair of weight, weight specification for 24 inch) 1090g



Handrims

Ellipse 3R, black

Item Code: QSF070212, QSS070212

Ellipse 3R handrims ensures pushing the wheelchair is easier. With its ergonomic curved shape and design, the Ellipse 3R mirrors the hand's natural shape for enhanced propelling. (Pair of weight, weight specification for 24 inch) 1270g



Handrims

Airy (Omobic), black

Item Code: QSF070105, QSS070105

Ultralight ergonomic handrim with optimized profile for different hand sizes. Anti-slip silicone insert ensures reliable grip—even when wet—and is easy to clean. 840g



Handrims

Flow (Omobic), Lizard grip, black

Item Code: QSF070112, QSS070112

Extremely lightweight ergonomic handrim with optimized profile for both small and large hands. Provides a secure grip—even when wet—and is available with multiple surface treatment options (including dipping) to meet different user needs, including tetraplegic users.
640g



Handrims

Handrim cover, rubber

Item Code: QSF070250, QSS070250

A rubber cover that is fitted over the handrim which can also be removed.
(Pair of weight, weight specification for 24 inch)
934g



Handrim mount, super narrow (7 mm)

Item Code: QSF070307, QSS070307

17 mm (070317) is the standard setting that most users will be used to. It is also possible to have the handrim mounted closer at 7 mm (070307) and spaced out at 27 mm (070327). (Is only an adjustment and therefore no weight)
0g



Handrim mount, narrow (17 mm)

Item Code: QSF070317, QSS070317

17 mm (070317) is the standard setting that most users will be used to. It is also possible to have the handrim mounted closer at 7 mm (070307) and spaced out at 27 mm (070327). (Is only an adjustment and therefore no weight)
0g



Handrim mount, wide (27 mm)

Item Code: QSF070327, QSS070327


17 mm (070317) is the standard setting that most users will be used to. It is also possible to have the handrim mounted closer at 7 mm (070307) and spaced out at 27 mm (070327). (Is only an adjustment and therefore no weight)
0g



0°

Item Code: QSF070600, QSS070600


Camber is the angle of rear wheels in relation to the ground, typically at 90° (0 camber). Negative camber (e.g., in sports wheelchairs) angles the wheels inward, reducing force needed to maneuver and improving lateral stability due to increased wheel spacing. However, it may be harder to navigate narrow spaces. (Weight included in basic chair weight.)
0g



2°

Item Code: QSF070602, QSS070602


Camber is the angle of rear wheels in relation to the ground, typically at 90° (0 camber). Negative camber (e.g., in sports wheelchairs) angles the wheels inward, reducing force needed to maneuver and improving lateral stability due to increased wheel spacing. However, it may be harder to navigate narrow spaces. (Weight included in basic chair weight.)
0g



4°

Item Code: QSF070604, QSS070604


Camber is the angle of rear wheels in relation to the ground, typically at 90° (0 camber). Negative camber (e.g., in sports wheelchairs) angles the wheels inward, reducing force needed to maneuver and improving lateral stability due to increased wheel spacing. However, it may be harder to navigate narrow spaces. (Weight included in basic chair weight.)
0g



Princess

Item Code: QSF070823, QSS070823

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Basketball

Item Code: QSF070802, QSS070802

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Basketball 2

Item Code: QSF070803, QSS070803

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Chrome

Item Code: QSF070806, QSS070806

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Feather

Item Code: QSF070811, QSS070811

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Pool

Item Code: QSF070833, QSS070833

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Zebra

Item Code: QSF070822, QSS070822

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Flying Water

Item Code: QSF070812, QSS070812

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Galactica

Item Code: QSF070813, QSS070813

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Graffiti

Item Code: QSF070814, QSS070814

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Happy Lions

Item Code: QSF070816, QSS070816

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Hot Wheels

Item Code: QSF070817, QSS070817

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Jungle Fever

Item Code: QSF070818, QSS070818

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Oriental

Item Code: QSF070820, QSS070820

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Pool

Item Code: QSF070822, QSS070822

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Monkeys

Item Code: QSF070819, QSS070819

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Pirates

Item Code: QSF070821, QSS070821

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Silk

Item Code: QSF070827, QSS070827


Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Space

Item Code: QSF070828, QSS070828

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Superhero

Item Code: QSF070829, QSS070829

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Tiger

Item Code: QSF070830, QSS070830


Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Transparent

Item Code: QSF070831, QSS070831

Spokeguards cover the rear wheel and thus protect against pinching fingers when accidentally coming into contact with the spokes of the wheels. The spokeguards come in large variety of designs. (Pair of weight)
780g



Standard brake

Item Code: QSF060001, QSS060001

Push to lock, spring loaded for ease of use, when engaged the lever lies out of the way making side transfers easy. (Pair of weight)
498g



Knee-lever brake

Item Code: QSF060002, QSS060002

Push to lock, lower forces to engage the brake, flip-down lever to prevent unintended opening by swing-away hangers and make room for transfers. (Pair of weight)
558g



Brake lever extension, knee-lever brake, fold-down

Item Code: QSF060010, QSS060010

This extension handle makes operating the knee lever wheel lock easier for users with restricted movement. (Pair of weight)
100g



Compact brake EVO

Item Code: QSF060023, QSS060023

This wheel lock is lightweight and low profile. It's out of the way for transfers but requires hand-function for engaging/ disengaging. (Pair of weight)
238g



Compact brake lightweight EVO

Item Code: QSF060024

This wheel lock is lightweight and low profile. It's out of the way for transfers but requires hand-function for engaging/ disengaging. (Pair of weight)
190g



Brake, pull to lock, brake lever, long

Item Code: QSF060009, QSS060009

Pull-to-lock brake with extended lever for increased reach and easier engagement.
422g




Aluminium compact brake (Omobic), black

Item Code: QSF060031, QSS060031

This compact aluminium scissor brake features an ergonomic grip and is also suitable for reduced hand strength. Its low-profile design prevents protrusion, making it easier to lift the wheelchair into a vehicle.

376g



Tool kit

Item Code: QSF090000, QSS090000

Contains the tools required for the adjustments available on the wheelchair.

218,5g



High pressure pump (0 - 11 bar)

Item Code: QSF090024, QSS090024

Most of our pneumatic tyres are high pressure cycling tyres and can be pressurised up to 140 psi (please refer to the recommended amount on the tyre). The high pressure pump can fill the tyre to this level although some users may require assistance to do so.

84g



Air pump

Item Code: QSF090025, QSS090025

A simple hand pump that can keep the tyre inflated to a reasonable pressure but is very difficult to reach the true pressure for the chair. Useful for the user to keep with them to keep punctures inflated while out and about until they can be repaired.

1148g



Back pack

Item Code: QSF090030, QSS090030

Designed with special hoops to hang from the push handles or backposts of a wheelchair. Please keep your belongings secure.

541g



Frame reinforcement kit for add on bike

Item Code: QSF090014

Reinforcement components designed to increase frame stability and durability when using a front add-on bike.

220g



Mobile phone pocket

Item Code: QSF090026, QSS090026

Holds the mobile at the front of the chair.

52g





Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

For further information on the full specification, options, and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.



Rev. A 04.2026



www.SunriseMedical.eu