

# OPTIMUS 2/ OPTIMUS 2 RS 2.322

## USER PROFILE

✓ For people using a manual wheelchair or walking aid indoors and who require more support outdoors

## HIGH-PERFORMANCE OUTDOOR POWER WHEELCHAIR (WITH FRONT-WHEEL DRIVE)

More powerful batteries and a top speed of up to 15 km/h make the all-terrain power wheelchair an experience.

- Ideal for long distances out of doors
- Optimised chassis geometry enables best-possible handling
- Single-wheel suspension and high climbing ability are ideal for high kerbs, higher slopes and challenging and rough ground
- Direct steering combined with a travel speed of 6, 10 or 15 km/h for high-level safety and reliability outdoors
- Larger batteries available for long range and outstanding climbing ability
- Extensive range of seats for optimally adapted fitting



RECARO seating systems available



Obstacle climbing ability up to 110 mm



Electric seat and back angle adjustment also for more complex fittings (optional)

R-Net electronics as standard

Puncture-proof tyres with studded tread for optimum traction and off-road capability

Drum brakes as standard



User weight up to 150 kg



Very powerful motor 950 W



Very high max. range between 50 and 77 km (depending on battery size)



# NEW

# DETAILS

OUTDOORS

OPTIMUS 2 2.322

### New 15 km/h software

Improved drive stability and range  
(up to 77 km in the 15 km/h version with 125 AH (20h) battery), approved for 18% gradient



### LED lighting

Bright front and rear lighting and new white front safety reflector

### Integrated legrest


Detachable, robust legrests



Mechanical height adjustment



Electric height adjustment



OPTIMUS 2  
MORE FREEDOM WITH  
THE RANGE MIRACLE

PICTURED:  
OPTIMUS 2  
PAGE 34