



OCTANE

SUB4

TITANIUM
TECHNOLOGY:
A NEW LEVEL OF DURABILITY



*Excluding wheels and cushion

OCTANE **SUB4**

THE WORLDS FIRST MADE TO MEASURE
TITANIUM WHEELCHAIR GUARANTEED
UNDER 4KG*

The Octane Sub4 elevates the daily wheelchair to a whole new level of performance and engineering down to the finest detail.

We took a classic design and redefined it even further by combining Sub4 technology with a hand-crafted made to measure titanium grade 9 frame. Resulting, in an ultra-lightweight and rigid wheelchair that is perfectly designed to fit you and your lifestyle.

TITANIUM GRADE 9:

- STRENGTH TO WEIGHT RATIO
- SMOOTH RIDING
- NO CORROSION
- DURABILITY





Frame finish options: polished or brushed



Sub4 technology – guaranteed under 4 kg*

KEY FEATURES

- Made to measure
- Less than 4 kg guaranteed (excluding wheels and cushion)
- Rigid and lightweight
- Titanium grade 9 frame
- Low profile castor arm design
- Airtech Sub4 technology
- Low profile cross braces
- Optimised single tube design
- Max user weight 110 kg



Low profile castor arm design



Single post axle

OCTANE **SUB4**

TECHNICAL DATA

OCTANE SUB4 // Titanium Wheelchair

Max. user weight	110 kg	Total length	880 mm
Frame version	Fully welded open frame	Total height	900 mm
Frame inset	Custom-made	Max. pivot width	880 mm
Frame angle	Custom-made	Max. safe slope	7° (for use of parking brakes)
Seat width	250 - 440 mm	Front wheel	3" / 4" / 5"
Seat depth	250 - 440 mm	Rear wheel	24" / 25"
Seat height front	380 - 520 mm	Wheelchair weight	From 6.2 kg
Seat height rear	360 - 500 mm	Dismanteled weight	From 3.6 kg
Backrest height	100 - 400 mm	Mass of heaviest part	2.1 kg
Backrest angle	Custom-made	Removable parts	Rear wheels (for transport only)
Centre of gravity (COG)	25 - 125 mm	Colours	- Frame finish options: Brushed frame or polished frame - 5 accent colours (decal on frame), 6 backrest upholstery stitching colours, 9 Spinerger colour options
Wheel camber	1° / 2° / 3°	Intended use as seat in motor vehicle	Not released as seat in motor vehicles
Total width	Max. 680 mm (SW + 240 mm)	Transport without user	Fix the wheelchair safe in the vehicle (e.g. car or plane)
		Intended use and environment	For people with limited mobility (children & adults); self-propelling or being pushed by an assistant, indoor and outdoor use.



Sunrise Medical ApS

Mårkærvej 5-9
2630 Taastrup
Tlf.: 70 22 43 49
E-mail: info@sunrisemedical.dk
www.SunriseMedical.dk

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.

