

Page 1 of 68 Report Reference No.: 713186041

TEST REPORT EN 12184:2014

TUV SUD Test report for

Electrically powered wheelchairs, scooters and their chargers -Requirements and test methods

	•		
Report reference No:	713186041		
Date of issue:	2020-12-15		
Project handler:	Matthias Müller		
Testing laboratory:	TÜV SÜD Product Service GmbH		
Address:	Masurenweg 1-3, 30163 Hanover, Germany		
Testing location:	as above		
Client:	Sunrise Medical GmbH		
Client number:	36070		
Address:	Kahlbachring 2-4, 69254 Malsch / Heidelberg, Germany		
Contact person	Mr. Gerhard Weis		
Standard::	This TUV SUD test report form is based on the following requirements: EN 12184:2014		
TRF originated by:	TUV SUD Product Service GmbH, Mr. Michael Kese, DiplIng. (product specialist)		
Copyright blank test report::	This test report is based on the content of the standard (see above). The test report considered selected clauses of the a.m. standard(s) and experience gained with product testing. It was prepared by TUV SUD Product Service GmbH.		
	TUV SUD Group takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.		
Scheme:	☐ TÜV Mark, ⊠ without certification		
Non-standard test method:	⊠ No □ Yes, see details under Summary		
National deviations:	N/A		
Number of pages (Report):	68		
Number of pages (Attachments):	N/A		
Compiled by: Matthias Müller	Approved by: Torsten Zimmer		
(+ signature)	(+ signature)		
SICN ID 4267	SIGN-ID 340341		
SIGN-ID 43670	16 12 2020		

15.12.2020



Page 2 of 68 Report Reference No.: 713186041

Test sample:	HAN-502600-1			
Type of test object:	Electrically driven wheelchair			
Trademark::	SUNRISE MEDICAL.			
Model and/or type reference::	Quickie Q300M Mini			
Rating(s):	See copy of marking plate			
Manufacturer:	Sunrise Medical GmbH			
Manufacturer number:	36070			
Address:	Kahlbachring 2-4, 69254 Malsch / Heidelberg, Germany			
Sub-contractors/ tests (clause):	N/A			
Name:	N/A			
Order description:				
	☐ Partial test according to manufacturer's specifications			
	☐ Preliminary test			
	☐ Spot check			
Date of order:	2020-04-06			
Date of receipt of test item::	: 2020-07-01			
Date(s) of performance of test:	: 2020-07-03 to 2020-12-15			
Test item particulars:				
Electrically powered wheelchair, type Quickie Q300M Mini				
Attachments: N/A				
General remarks:				
"(see remark #)" refers to a remark appended to the report. "(see appended table)" refers to a table appended to the report. Throughout this report a comma is used as the decimal separator. The test results presented in this report relate only to the object tested.				
This report shall not be reproduced except in full without the written approval of the testing laboratory.				



Page 3 of 68 Report Reference No.: 713186041

Summary of testing:			
The product fulfils the requirments	of EN 12184:2014		
Additional information on Non-sta	ndard test method(s)		
Sub cl	::		
Page	::		
Rational	::		
☐ deviation(s) found / product is ☐ no deviations found / product If additional information is necessary.	is in compliance with th	·	
Copy of marking plate:	Sunrise Medica Kahlbachring 2	2-4 2020-12-08	3
Type: POL	69254 Malsch/l	0 A 116-20-415-34 (1-10-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-	_
100 - 100 -			- :
Q300M MINI	ISO 7176-19:2008	FIN-Nr: 3M20128328	
136 kg 285 kg	10 km/h 8°	60/250/100kg	